

Hate your crutches? Limited by your knee scooter? THE NEXT GENERATION OF HANDS-FREE MOBILITY IS HERE!

Easily manage you daily activities, play with your kids, hold your baby, get up and down stairs, use public transit, work, exercise, or go shopping!

:: NEW & IMPROVED



- Lateral offset, catamaran rocker tread, selfstabilizing foot makes balancing easy and intuitive. Rotates 180 for right or left leguse.
- Proprietary glass reinforced thermoplastic, paired with engineered forms, keeps the device light and sturdy.



- Comfortable woven straps have fabric covered foam pads, easy infinite adjustment and quick release buckles.
- Donning and doffing takes seconds.



Integral handle makes learning easier and allows short trips without strapping in. Triangulated alloy tube frame enhances user stability.



- Ergonomic knee platform shape provides greater comfort and control.
- Independent upper and lower leg height adjustments makes fitment fast and easy.



₩ TESTIMONIALS

family, and my sanity despite a nasty injury.

"As a General Surgeon, I spend my working hours standing at an operating table,

impossible with crutches...and have been on my iWALKFREE...for nearly six months. This device is one of the most fantastic, original, insightful and important

"This product has changed my outlook. I feel better now that I can work and be more useful at home with simple things that I could not do on crutches."

"[iWALK] provides many benefits for the patient...early mobilization and help[ing] to decrease loss of both bone density and muscle strength.

medical innovations of recent years and has allowed me to maintain my job, my

- Rigid, independently adjustable thigh cradles provide enhanced stability and simple tool-free size adjustments.

-Richard W. Shallman MD. WA

-Janet Cleave H.B.Sc., B.Sc., P.T.

ELIMINATE HAND / UNDER ARM / WRIST PAIN PEASIER THAN CRUTCHES The state of the s **:** LIFESTYLE

GET YOUR LIFE BACK

EASY TO LEARN

Lightweight, durable, comfortable and inexpensive, the world's best selling hands-free crutch just got even better.

Fits heights 4'10" - 6'6", 275 lbs (147-198 cm, 125 kg)

For lower leg injuries only. Should not be used





EASY TO USE

AFFORDABLE















COMFORTABLE PADDING

ADJUSTABLE THIGH SUPPORTS (18" - 27") (46-69cm)

EASY ADJUST FABRIC STRAPS

QUICK ON/OFF BUCKLES

REPLACEABLE TREADS




Fits heights 4'10" – 6'6", 275 lbs (147-198 cm, 125 kg)

INSTANT

PUSH BUTTON

HEIGHT ADJUST

FEATURES -





if balance or leg strength is compromised.





For lower leg injuries including:

Achilles injuries, foot / ankle sprains and fractures, foot ulcers, tibia / fibula fractures Most lower leg non-weight bearing injuries

BENEFITS:

- Improves mobility during rehabilitation
- Reduce muscle atrophy
- Reduces post-recovery rehabiliation time
- Affordable
- Covered by many insurance plans
- Eliminates pain and fatigue of conventional crutches
- Boosts user independence, morale and confidence
- Ability to perform common tasks that are impossible with standard crutches & knee scooters
- Easier to learn than conventional crutches
- Fully adjustable Fits users 4'10" to 6'6"
- Adapts to left or right leg use in seconds without tools



Easy as 1-2-3



Easy, quick, tool free



Self fit in minutes



One small step for man
One giant leap for mobility



iwalk-free.com

The world's original hands free crutch Universes

Distributor information here

iWALKFree, Inc. 956513 C7 Mulmur, Mansfield, ON LON1MO, Canada 877 268-1466 sales@iwalk-free.com

LEARN MORE HERE!

For lower leg non-weight bearing injuries

GET YOUR LIFE BACK



HANDS FREE PAIN FREE KNEE CRUTCH



For more information go to iwalk-free.com









