



135
YEARS
EST. 1888

Celebrating
135 Years of
Commitment
to Respiratory
Care

drive
 **DeVilbiss**
HEALTHCARE

135 Years of Respiratory Excellence



Maximize Patient Outcomes

With the Most Innovative, Dependable Products on the Market

From pioneering early respiratory technology to aiding in the pandemic, Drive DeVilbiss Healthcare has achieved significant milestones in the respiratory field across its 135-year history.

With more than 135 years of commitment to respiratory care and a full catalog of respiratory products, we've weathered countless changes to continue offering our customers the products they need. We look forward to staying by your side for the road ahead.



iGo[®]2 Portable Oxygen Concentrator With Optional Bluetooth Connectivity

- Item # 125D iGo[®]2 Portable Oxygen Concentrator
- Item # 125D-XB iGo[®]2 Portable Oxygen Concentrator with Extra Battery Pack

- Settings 1-5 PulseDose[®], 6 SmartDose[®]
- SmartDose[®] Auto-Adjusting Technology included
- Rugged design, including protective overmold
- Weighs less than 5 lbs.
- Most sensitive conserver trigger on the POC market
- 5-year warranty (3-year compressor, 1-year sieve beds & battery)



Proudly
Assembled in
The USA

[Learn More](#)

- Item # 125D-BT iGo[®]2 Portable Oxygen Concentrator with Bluetooth
- Item # 125D-BT-XB iGo[®]2 Portable Oxygen Concentrator with Bluetooth and Extra Battery Pack

Same features as base model with:

- Bluetooth Connectivity
- SmartLinkO2[®] App
- Remote Fleet Monitoring
- SmartLinkO2[®] Cloud



Proudly
Assembled in
The USA

[Learn More](#)

135 Years of Respiratory Excellence



5L & 10L Concentrator

Item # 525DS 5L Concentrator

- Built upon a legacy of reliable, concentrator models, the 5 Liter Oxygen Concentrators are designed with enhanced quality, durability, performance and simplicity in mind
- Exclusive DeVilbiss OSD® (Oxygen Sensing Device), standard on every unit, ensures patient safety and reliability for longer service intervals
- Readily accessible patient controls, protected cannula fitting and recessed humidifier nook to prevent damage
- Alarms: Audible and visual high/low pressure, low flow, low oxygen, power fail, Oxygen Sensing Device

[Learn More](#)



Proudly
Assembled
in the USA

Item # 1025DS 10L Concentrator

- Proudly assembled in the USA with the same shell as the 5 Liter, the 10 Liter Concentrator is also compatible with cylinder transfill systems
- With a weight of 42 lbs. and an identical footprint to the 525 unit, the 1025 is substantially smaller and lighter than competitive 10 Liter concentrators
- Capable of delivering 87% to 96% of O₂ purity from 2 to 10 LPM, the 10 Liter Concentrator is an ideal investment for your patients and business alike
- Exclusive DeVilbiss OSD® (Oxygen Sensing Device), standard on every unit, monitors the oxygen produced and ensures patient safety and reliability for longer service intervals

[Learn More](#)



Proudly
Assembled
in the USA

Bonsai® Velocity Pneumatic Oxygen Conserver

Item # OM-812

- Delivers a uniform oxygen pulse at increments: 1, 2, 3, 4, 5 and 6
- Two continuous flow (CF) settings preset @ 2 LPM and 4 LPM
- Delivers up to a 6:1 savings ratio from 14 – 40 breaths per minute
- Operating pressure range: 500-3000 PSI
- Light weight and quiet operation
- Easy-to-operate single selector knob



Proudly
Assembled
in the USA

[Learn More](#)



Evolution™ Electronic Oxygen Conserver

Item # OM-900

- Delivers a uniform oxygen pulse at increments: 1, 2, 3, 4, 5, 6, 7 Plus CF preset @ 2 LPM
- Minimum 2-year battery life* utilizing 2 AA batteries
- Delivers a 5:1 savings ratio from 14-40 breaths per minute
- Senses a breath/delivers oxygen in the first half of the inspiratory cycle
- Operating pressure range: 500-3000 PSI
- Weighs only 14.9 ounces with batteries
- Single push button selects the equivalent LPM dosage settings
- Auto shut-off feature conserves battery life



Proudly
Assembled
in the USA

[Learn More](#)



SmartDose® Mini Auto-Adjusting Oxygen Converter

Item # CTOX-MN02

- Delivers a uniform oxygen pulse at increments: 1, 2, 3, 4, 5 plus CF present @ 2 LPM
- Continuously monitors breathing patterns and adjusts oxygen dose to match activity level
- Achieves maximum conservation efficiency by automatically adjusting the dose down when less oxygen is needed
- Allows delivered dose to increase up to 2 settings during patient exertion up to a maximum pulse dose setting of 7
- Operating pressure range: 500 – 3000 PSI
- Weighs just 15.5 ounces with batteries
- Auto shut-off feature conserves battery life



Proudly
Assembled
in the USA

[Learn More](#)



iFill® Personal Oxygen Station

Item # 535D

- Offers a clinically focused and operationally efficient non-delivery oxygen solution
- Makes it easy and convenient for your customers to fill portable oxygen cylinders in their homes and gives them greater freedom to come and go as they please
- Fills proprietary iFill cylinders - Standard 870 Post Valve or Integrated Continuous Flow regulator
- Easily tops off partially filled cylinders and automatically shuts off when filling is complete
- Control panel light signals when cylinder is filling and when it's full



Proudly
Assembled
in the USA

[Learn More](#)



Pulmo-Aide® Compact Compressor Nebulizer

Item # 3655D With Disposable Neb Kit

- The DeVilbiss Pulmo-Aide® Compact Compressor Nebulizer System was built to last with high-quality metal mechanical parts, ensuring quality therapy for an extended period of time
- Small footprint
- Small, compact design
- Convenient carrying handle

[Learn More](#)



Proudly
Assembled
in the USA



PulmoMate® Compressor Nebulizer

Item # 4650D With Disposable Neb Kit

- The DeVilbiss PulmoMate® features an updated compressor/motor for long-term durability and performance
- Consistent performance over the life of the unit
- Metal frame motor construction
- Integrated handle for easy carrying
- Large built-in compartment for storing nebulizers and tubing
- Intuitive design is easy-to-use and maintain

[Learn More](#)



Proudly
Assembled
in the USA



Pulmo-Aide® Compressor Nebulizer

Item # 5650D With Disposable Neb Kit

- The DeVilbiss Pulmo-Aide® Compressor Nebulizer System is the only approved compressor for use with TOBI (Tobramycin)
- One of the longest-lasting compressors on the market with an average product life of 10 years
- Unmatched quality and durability – great for nursing homes and long-term users
- Metal frame motor construction
- Large, built-in compartment for storing nebulizers and tubing

[Learn More](#)



Proudly
Assembled
in the USA



Power Neb Ultra

Item # 18080 With Disposable Neb Kit

Item # 18081 With Disposable & Reusable Neb Kit

- A competitive price point for providers as one of the smallest compressor nebulizers at just 3.2 lbs
- Built-in neb kit holder for hands-free use
- Powerful piston pump is efficient, dependable and quiet
- Delivers consistent particle size of <5 microns and fast treatment times

[Learn More](#)



Pacifica Elite

Item # 18070 With Disposable Neb Kit

Item # 18071 With Disposable & Reusable Neb Kit

- This retail-packaged compressor neb is a workhorse with higher liter flows, shorter treatment times and quieter operation than similar economy compressor nebs.
- Ultra-compact design makes it extremely portable
- Powerful piston pump is efficient, dependable and quiet
- Delivers consistent particle size of <5 microns and fast treatment times

[Learn More](#)



Retail
Packaged



Vacu-Aide® Suction Machine

Item # 7325P-D With Battery, Internal Filter & Carry Bag

Item # 7325P-D-EXF With Battery, External Filter & Carry Bag

- Compact and lightweight
- Improved, longer-lasting lithium ion battery
- Best in class warranty: **two-year**, excluding container, **one-year** lithium ion battery
- Quiet, <59 dBA
- Standard canister, tubing and accessories
- Carry bag included with 7325P series

[Learn More](#)



Proudly
Assembled
in the USA





CORPORATE OFFICE  99 Seaview Blvd., Port Washington, NY 11050  877-374-8763  www.drivemedical.com

CANADA OFFICE  375 New Huntington Rd., Unit 1, Vaughan, Ontario L4H 5B5  866-752-1363

©2023 Medical Depot, Inc. dba Drive DeVilbiss Healthcare, 99 Seaview Boulevard, Port Washington, NY 11050. All trademarks used in association with the sale of products of Drive DeVilbiss are trademarks owned by Medical Depot, Inc. All other trademarks, trade names, service marks and logos referenced herein belong to their respective companies. Refer to product instruction guide for warranty details.